

HERizon™

Detect, Diagnose, and Treat Uterine Cavity Conditions



HERizon Hysto-Kit: All-in-one, single-use hysteroscopy kit

Provides essential supplies to perform diagnostic and therapeutic hysteroscopy procedures to help reduce setup time and aid with inventory management. Each sterilized kit includes:

- Back table cover
- Towel
- Spinal needle
- Irrigation tubing
- Alcohol prep pad
- 18g needle
- Prep and swab stick
- 10cc syringe
- Sterile glove
- Under buttock drape with collection pouch
- Wooden cotton-tip applicator
- Suction tubing
- Extension tube set

No cervical dilation


Minimal anesthesia

Reduce cross contamination

Compact office foot print

minerva[®]
SURGICAL

Hysteroscopy Accessories - Scissors, Forceps and Graspers



Scissors



Alligator



Foreign Body Grasper



Biopsy

To order contact Minerva Surgical Customer Service

By phone: call **855.646.7874**, and press for option 1, or email: customerservice@minervasurgical.com, or fax 866.465.2875

Product Name	Catalog Number	Quantity	Description
Medical Endoscope Image Processor	JY-MIP-1000	1 Unit	Camera and light source control unit
Disposable Electronic Hysteroscope	JY-DEH-001	1 Unit	Disposable 4.8 mm hysteroscope w/operative port, orders sold in multiples of 10 units
Disposable Endoscopic Forceps - Foreign Body Grasper	JY-C-350-18-S	1 Unit	For grasping a foreign body and removal under hysteroscopic guidance, orders sold in multiples of 10 units
Disposable Endoscopic Forceps - Alligator	JY-CC-350-18-S	1 Unit	For grasping, manipulating, and retrieving tissue, orders sold in multiples of 10 units
Disposable Endoscopic Forceps - Biopsy	JY-PK-350-18-S	1 Unit	For tissue biopsy or sampling under hysteroscopic guidance, orders sold in multiples of 10 units
Disposable Endoscopic Scissors	JY-JD-17-350-B	1 Unit	For tissue cutting under hysteroscopic guidance, orders sold in multiples of 10 units
HERizon Hysto-Kit	HERHYSTOKIT-1	1 Unit	Single-use hysteroscopy kit (e.g., drapes, tubing, towels, syringe, etc.), orders must be multiples of 5 units

Intended use: The Medical Endoscope Image Processing System is intended to be used for viewing of the adult cervical canal and uterine cavity for the purpose of performing diagnostic and operative procedures.

The Disposable Endoscopic Forceps is intended for grasping tissue and foreign bodies during endoscopic procedures.

The Disposable Endoscopic Scissors are designed to cut tissue during endoscopic procedures.

Contraindications: This product is not suitable for people with acute pelvic inflammatory disease. Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent: inability to distend the uterus, cervical stenosis, cervical/vaginal infection, uterine bleeding or menses, known pregnancy, invasive carcinoma of the cervix, recent uterine perforation, medical contraindication or intolerance to anesthesia.

Warnings: Failure to follow any instructions or to heed any Warnings or Precautions could result in serious patient injury. The system should only be used by physicians trained in hysteroscopy and hysteroscopic surgery using powered instruments. Healthy tissue can be injured, e.g., perforation. Use every available means to avoid such injury.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. The physician using the system must be trained in hysteroscopy. Prior to use, consult the Instructions for Use and/or Operators Manual. For Indications and a complete list of contraindications, potential adverse events, warnings and precautions.



Product manufactured and packaged
for Minerva Surgical, Inc.
By Jiangsu Jiyuan Medical Technology Co., Ltd.
4255 Burton Drive, Santa Clara, CA 95054 USA

www.minervasurgical.com
© 2025 Minerva Surgical and its affiliates.
All rights reserved

minerva
SURGICAL
K0309 Rev. C